Reprint from
Selected Topics in Clinical Enzymology, Vol. 2
Editors: David M. Goldberg, Mario Werner
© 1984 Walter de Gruyter & Co., Berlin · New York-Printed in Germany

DETERMINATION OF PSEUDOCHOLINESTERASE ACTIVITY BY AN ENZYME COUPLING METHOD IN THE SERUM OF PATIENTS WITH LIVER DISEASE

Kwang-Jen Hsiao, Shew-Jen Wu, Yasushi Kasahara, Wen-Kuo Ting

Clinical Biochemistry Research Laboratory, Department of Medical Research and Department of Laboratory, Veterans General Hospital, Taipei, Taiwan 112; and Division of Clinical Chemistry, Research Laboratories, Fujirebio Inc., Komiya, Hachioji, Tokyo 192, Japan.

Abstract

The serum pseudocholinesterase (EC 3.1.1.8, acylcholine acylhydrolase) activity was determined at 37°C and pH 8.2 by using a stable substrate p-hydroxybenzoylcholine as the initiator and p-hydroxybenzoate hydroxylase (EC 1.14.13.2) as the The reaction was measured by continuously coupling enzyme. monitoring the decrease of NADPH at 340nm. The within-run and between-run precisions of the test were 0.4-2.6% (CV) and 1.2-1.9%, respectively. The test was linear at least up to 600U/L and had the sensitivity of 0.001\DeltaA/min at 9.7U/L. serum pseudocholinesterase reference range was estimated to be 156-418U/L (95% range; mean = 251U/L) in 200 healthy There was no significant difference between male and adults. The serum pseudocholinesterase activity of patients with liver diseases was significantly lower than that of normal controls, and in 90-95% of the cirrhosis and hepatoma patients was below the lower normal limit. This method for determination of serum pseudocholinesterase activity is suitable as a diagnostic aid for liver diseases.

Introduction

Pseudocholinesterase (EC 3.1.1.8, acylcholine acylhydrolase) is an enzyme which catalyzes the hydrolysis of various chol-